



## Fit for a Genius Discover Einstein™

Meet the Einstein 3D printer series. Meticulously designed for dental professionals, Einstein delivers the key features in 3D printing: Accuracy, Speed and Versatility.

---

## The Secret to Einstein is in its DNA

---

**Proprietary NanoFit 385** technology ensures that every application features stunning clarity, perfect fit and an impressively natural looking finish.

**Equipped with HyperPrint,** harnessing the power of heat and a closed-loop process, fabricate a variety of dental applications with ease and up to 50% faster speed than its predecessor.

**Versatility** in materials, from dentures to models and everything in between, designed to maximize productivity for your lab or office.

---

**Transform your workflow from fit to finish with Einstein**

---



Properties	Einstein
HyperPrint	✓
<b>NanoFit 385</b>	✓
Technology	Digital Light Processing (DLP)
Wavelength	385nm
Build Envelope	7.48" x 4.21" x 4"
XY Native Resolution	99µm
Enhanced XY Resolution (Pixel Tuning)	65µm
Dynamic Layer Thickness (Z Resolution)	25-50-75-100-150µm
Electrical Requirements	100–250V AC, 5A, 50Hz/60Hz
Footprint	16" x 16" x 27.5"
Weight	75lbs
Included Warranty	1-year
Printer SKU#	MIC-55-1000
Technical Support	Monday through Friday, 8am–8pm Eastern