Desktop Health



Introducing

Model X & Model Z



Model X with its 50-micron build layer allows for exceptionally detailed crowns, bridges with removable dies, implant cases, and diagnostic and orthodontic models.

Flexural Strength: 90 Mpa

Flexural Modulus: 2.2 GPa

Elongation at Break: 5.5%











High Speed Meet Premium Precision



Model Z is as effective as it is fast. This low viscosity model material is ideal for orthodontic, diagnostic, and crown and bridge models requiring smooth surface resolution, easy clean up, and a superior fit and finish.

Flexural Strength: 90 Mpa

Flexural Modulus: 2.5 GPa

Elongation at Break: 10%





DesktopHealth.com @DesktopHealth (866) 344-3578